

Product Code . EDS-E-10592

Force on a Conductor Apparatus



Description

When the axle is placed on the rails between the poles of the magnet, and a power supply unit is connected, the axle is repelled and rolls along the rails away from the centre of the magnetic field.

Jainco Export comprising a strong U-shaped magnet and a pair of brass rails with 4mm socket terminals.

A brass axle with plastic discs is free to roll along the rails and completes the Electrical contact between them.

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha, Haryana
Direct Contact Details ☎ +91-98173-19615 ✉ sales@educational-equipments.com
💻 www.educational-equipments.com