



Product Code . EDS-H-10677

Conductivity Apparatus

Description

Jainco Export Rates of melting compared by placing wax on the ends of the different rods thus demonstrates rates of conductivity in the different metals.

Jainco Export Five rods, one each of aluminium, copper, brass, zinc and iron, embedded along one side of a metal tank, size 150 x 90 x 100mm.

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha, Haryana
Direct Contact Details ☎ +91-98173-19615 ✉ sales@educational-equipments.com
💻 www.educational-equipments.com