



Product Code . EDS-H-10686

Conductivity of Metal Apparatus

Description

Gentle heating at the centre melts wax placed in the hollows at the end of each strip.

Jainco Export strips of copper, brass, steel and aluminium mounted on a wooden ring and meeting in the centre.

Rates of conductivity in the different metals can be compared by the difference in time taken for each cup of wax to melt.

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha, Haryana
Direct Contact Details ☎ +91-98173-19615 ✉ sales@educational-equipments.com
💻 www.educational-equipments.com