

Educational Lab Equipments





Product Code . EDS-PLE-10814



Photosynthesis Apparatus

Description

The gas evolved can be measured using the calibrated capillary tube and drawn off for analysis using the plastic syringes.

Jainco Export for the measurement and collection of gas evolved by aquatic plants during photosynthesis.

Mounted on a plastic panel, 190mm x 150mm.

Educational Lab Equipments,

#449, HSIIDC, Industrial Area, Saha, Haryana

Direct Contact Details \(\subseteq +91-98173-19615 \) Sales@educational-equipments.com

www.educational-equipments.com