



Product Code . EDS-PLE-10816

Potometer Apparatus on White Stand

Description

When transpiration occurs, the bubble within the capillary moves.

A syringe is used to zero the scale.

Allows the students to measure transpiration with minimal risk of breakages.

Jainco Export used to show transpiration which occurs in small plant shoots.

The apparatus is submerged and the specimen placed in the tube.

The apparatus has only one glass component (the capillary tube) and is mounted on a strong backing plate.

The apparatus can also be used to measure humidity by connecting a porous pot.

Dimensions: 195x150x75mm

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha, Haryana
Direct Contact Details ☎ +91-98173-19615 ✉ sales@educational-equipments.com
💻 www.educational-equipments.com