

Jainco Export Horizontal support rod and bosshead included.

A soil sample is placed on the perforated plate and a 220V heating lamp is held above the soil, driving the organisms down through the funnel into a collecting jar.

Supplied with coarse and fine mesh.

Specification:

Sieve diameter 170mm

Sieve depth 65mm

Requires a low wattage bulb (Bayonet fitting 22mm diameter, 40 watt or less).

Funnel diameter 175mm

Educational Lab Equipments, #449, HSIIDC, Industrial Area, Saha, Haryana Direct Contact Details 4 +91-98173-19615 Sales@educational-equipments.com www.educational-equipments.com